

## **Certified Reference Material - Acidity, waste water**

Art. ID	ERA-915
Unit	each
Deliverydetails	No Dangerous Good /not restricted

### Description

One 250 mL whole-volume bottle is ready to analyze. Designed for use with titrimetric methods to a pH endpoint of 8.3. The certified values vary from batch to batch. Please ask for the certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Acidity as CaCO <sub>3</sub>		650 - 1800	mg/L		