

Certified Reference Material - Boron, waste water (WasteWatR TM)

Art. ID ERA-919
Unit each
Deliverydetails No Dangerous Good /not restricted

Description

One unpreserved 60 mL poly bottle yields in excess of 2 liters after dilution. Designed for colorimetric methods. The certified values vary from batch to batch. Please ask for the certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Boron (B)	[7440-42-8]	0,8 - 2	mg/L		