

Certified Reference Material - Uranium, potable water (PotableWatR TM)

Art. ID	ERA-930
Unit	each
Deliverydetails	ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN3264 III E1

Description

One 15 mL screw-cap vial yields up to 2 liters after dilution. Use with ICP-MS methods. The certified values vary from batch to batch. Please ask for the certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Uranium (U)	[7440-61-1]	3 - 104	µg/L		