

Certified Reference Material - Corrosivity, potable water (PotableWatR TM)

| | |
|-----------------|-----------------------------------|
| Art. ID | ERA-980 |
| Unit | each |
| Deliverydetails | No Dangerous Good /not restricted |

Description

One 500 mL whole-volume bottle is ready to analyze for corrosivity, calcium carbonate saturation and Langelier saturation index. The certified values vary from batch to batch. Please ask for the certificate.

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|------------|---------------------|--------------------|--------|--------|
| | Corrosivity | | -4 | to +4 SI u nits | | |