

## **Certified Reference Material - Corrosivity, potable water (PotableWatR TM)**

Art. ID	ERA-980
Unit	each
Deliverydetails	No Dangerous Good /not restricted

### Description

One 500 mL whole-volume bottle is ready to analyze for corrosivity, calcium carbonate saturation and Langelier saturation index. The certified values vary from batch to batch. Please ask for the certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Corrosivity		-4	to +4 SI u nits		