

## **Certified Reference Material - Cyanide, potable water (PotableWatR TM)**

|                 |   |
|-----------------|---|
| Art. ID         | ERA-983   |
| Unit            | 15 mL   |
| Deliverydetails | ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1824 III E1 |

### Description

One 15 mL screw-cap vial yields up to 2 liters after dilution. The certified values vary from batch to batch. Please ask for the certificate.

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|------------|---------------------|------|--------|--------|
|                  | Cyanide           |            | 0,1 - 0,5           | mg/L |        |        |
|                  | Free cyanide      |            | 0,1 - 0,5           | mg/L |        |        |
|                  | Total cyanide     |            | 0,1 - 0,5           | mg/L |        |        |