

## **Certified Reference Material - Cyanide, potable water (PotableWatR TM)**

Art. ID	ERA-983
Unit	15 mL
Deliverydetails	ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1824 III E1

### Description

One 15 mL screw-cap vial yields up to 2 liters after dilution. The certified values vary from batch to batch. Please ask for the certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Cyanide		0,1 - 0,5	mg/L		
	Free cyanide		0,1 - 0,5	mg/L		
	Total cyanide		0,1 - 0,5	mg/L		