

Proficiency Testing - Organochlorine Pesticides on Polyurethane foam

Art. ID ERA-1011
Unit each
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

One 2 mL flame-sealed ampule plus one polyurethane foam (PUF). Contains at least 16 analytes, randomly selected from the list below, at 0.5-20 µg/sample after preparation.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	4,4'-DDD	[72-54-8]	0,5 - 20	µg/sample		
	4,4'-DDE	[72-55-9]	0,5 - 20	µg/sample		
	4,4'-DDT	[50-29-3]	0,5 - 20	µg/sample		
	Aldrin	[309-00-2]	0,5 - 20	µg/sample		
	Dieldrin	[60-57-1]	0,5 - 20	µg/sample		
	Endosulfan-I (alpha)	[959-98-8]	0,5 - 20	µg/sample		
	Endosulfan-II (beta)	[33213-65-9]	0,5 - 20	µg/sample		
	Endosulfan-sulfate	[1031-07-8]	0,5 - 20	µg/sample		
	Endrin	[72-20-8]	0,5 - 20	µg/sample		
	Endrin aldehyde	[7421-93-4]	0,5 - 20	µg/sample		
	Endrin ketone	[53494-70-5]	0,5 - 20	µg/sample		
	Heptachlor	[76-44-8]	0,5 - 20	µg/sample		
	Heptachlor Epoxide (beta)	[1024-57-3]	0,5 - 20	µg/sample		
	Methoxychlor	[72-43-5]	0,5 - 20	µg/sample		
	alpha-HCH	[319-84-6]	0,5 - 20	µg/sample		
	cis-Chlordane (alpha)	[5103-71-9]	0,5 - 20	µg/sample		
	beta-HCH	[319-85-7]	0,5 - 20	µg/sample		
	delta-HCH	[319-86-8]	0,5 - 20	µg/sample		
	gamma-HCH (Lindane)	[58-89-9]	0,5 - 20	µg/sample		
	trans-Chlordane (gamma)	[5103-74-2]	0,5 - 20	µg/sample		