

## **Proficiency Testing - PCBs on Polyurethane foam**

Art. ID ERA-1012  
Unit each  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

### Description

One 2 mL flame-sealed ampule plus one polyurethane foam (PUF). Contains one Aroclor, randomly selected from the list below, at 2-10 µg/sample after preparation.

| Text/Information | Analyte/Parameter | CAS number   | Concentration/Value | Unit      | Method | Source |
|------------------|-------------------|--------------|---------------------|-----------|--------|--------|
|                  | Aroclor 1016      | [12674-11-2] | 2 - 10              | µg/sample |        |        |
|                  | Aroclor 1221      | [11104-28-2] | 2 - 10              | µg/sample |        |        |
|                  | Aroclor 1232      | [11141-16-5] | 2 - 10              | µg/sample |        |        |
|                  | Aroclor 1242      | [53469-21-9] | 2 - 10              | µg/sample |        |        |
|                  | Aroclor 1248      | [12672-29-6] | 2 - 10              | µg/sample |        |        |
|                  | Aroclor 1254      | [11097-69-1] | 2 - 10              | µg/sample |        |        |
|                  | Aroclor 1260      | [11096-82-5] | 2 - 10              | µg/sample |        |        |