

## **Proficiency Testing - PCBs on Polyurethane foam**

Art. ID ERA-1012  
Unit each  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

### Description

One 2 mL flame-sealed ampule plus one polyurethane foam (PUF). Contains one Aroclor, randomly selected from the list below, at 2-10 µg/sample after preparation.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aroclor 1016	[12674-11-2]	2 - 10	µg/sample		
	Aroclor 1221	[11104-28-2]	2 - 10	µg/sample		
	Aroclor 1232	[11141-16-5]	2 - 10	µg/sample		
	Aroclor 1242	[53469-21-9]	2 - 10	µg/sample		
	Aroclor 1248	[12672-29-6]	2 - 10	µg/sample		
	Aroclor 1254	[11097-69-1]	2 - 10	µg/sample		
	Aroclor 1260	[11096-82-5]	2 - 10	µg/sample		