

Proficiency Testing - Metals on Filter Paper

Art. ID	ERA-1025
Unit	each
Deliverydetails	No Dangerous Good /not restricted

Description

One filter paper sample is ready to be analyzed.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Antimony (Sb)	[7440-36-0]	30 - 1200	µg/filter		
	Arsenic (As)	[7440-38-2]	30 - 1200	µg/filter		
	Barium (Ba)	[7440-39-3]	30 - 1200	µg/filter		
	Beryllium (Be)	[7440-41-7]	30 - 1200	µg/filter		
	Cadmium (Cd)	[7440-43-9]	30 - 1200	µg/filter		
	Chromium (Cr)	[7440-47-3]	30 - 1200	µg/filter		
	Cobalt (Co)	[7440-48-4]	30 - 1200	µg/filter		
	Copper (Cu)	[7440-50-8]	30 - 1200	µg/filter		
	Lead (Pb)	[7439-92-1]	30 - 1200	µg/filter		
	Manganese (Mn)	[7439-96-5]	30 - 1200	µg/filter		
	Nickel (Ni)	[7440-02-0]	30 - 1200	µg/filter		
	Phosphorous		30 - 1200	µg/filter		
	Selenium (Se)	[7782-49-2]	30 - 1200	µg/filter		
	Silver (Ag)	[7440-22-4]	30 - 1200	µg/filter		
	Thallium (Tl)	[7440-28-0]	30 - 1200	µg/filter		
	Zinc (Zn)	[7440-66-6]	30 - 1200	µg/filter		