

## **Proficiency Testing - Solids, WP, waste water**

Art. ID	ERA-241
Unit	each
Deliverydetails	No Dangerous Good /not restricted

### Description

One 500 mL whole-volume bottle is ready to analyze. The CRM is also certified for pH. For a pH PT, order Cat. #577.  
Methods: EPA 160.1 EPA 160.2 EPA 160.3 SM 18/19thED 2540 C SM 18th ED 2540 B SM 18th ED 2540 D

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Total Suspended Solids	[n/a]	20 - 100	mg/L		
	Total dissolved solids at 180 °C		140 - 800	mg/L		
	Total solids at 105 °C		140 - 800	mg/L		