

## **Proficiency Testing - Solids, Concentrate, waste water**

Art. ID	ERA-4030
Unit	24 mL
Deliverydetails	No Dangerous Good /not restricted

### Description

One 23 mL screw-cap vial with a powder yields 1 liter after dilution.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Total suspended solids (TSS)		20 - 100	mg/L		
	Total dissolved solids at 180 °C		140 - 800	mg/L		
	Total solids at 105 °C		140 - 800	mg/L		