

## **Proficiency Testing - Texas Low-Level Fuels (TPH) in Water**

Art. ID                      ERA-476  
Unit                         each  
Deliverydetails            No Dangerous Good /not restricted

### Description

One 2 mL flame-sealed ampule yields in excess of 200 mL after dilution. Contains unleaded regular gasoline and No. 2 Diesel Fuel resulting in TPH in the range 5-10 mg/L.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Total petroleum hydrocarbons (TPH)		5,0 - 10,0	mg/L		