

Proficiency Testing - HEM / SGT-HEM, waste water

Art. ID ERA-489
Unit each
Deliverydetails No Dangerous Good /not restricted

Description

One 5 mL flame-sealed ampule yields up to 2 liters after dilution. Contains both hexadecane and stearic acid. Note: If a NELAC compliant PT is required, use Cat. #582 or Cat. #4120.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	HEM		5 - 100	mg/L		
	SGT-HEM		5 - 100	mg/L		