

## **Proficiency Testing - Solids, Concentrate, WS, potable water**

Art. ID                      ERA-5150  
Unit                         each  
Deliverydetails            No Dangerous Good /not restricted

### Description

One 23 mL screw-cap vial with a powder yields 1 liter after dilution.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Total suspended solids (TSS)		23 - 100	mg/L		
	Total dissolved solids at 180 °C		200 - 450	mg/L		
	Total solids at 105 °C		223 - 550	mg/L		