

## **Proficiency Testing - o-Phosphate Nutrients, potable water**

Art. ID ERA-558  
Unit each  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1789 III E1

### Description

One 15 mL screw-cap vial yields up to 2 liters after dilution.

| Text/Information | Analyte/Parameter    | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|----------------------|------------|---------------------|------|--------|--------|
|                  | Ortho-Phosphate as P |            | 0,5 - 5,5           | mg/L |        |        |