

## **Proficiency Testing - Tin & Titanium, waste water**

Art. ID	ERA-573
Unit	each
Deliverydetails	No Dangerous Good /not restricted

### Description

One 15 mL screw-cap vial yields up to 1 liter after dilution. Use with AA, ICP-OES or ICP-MS methods.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Tin (Sn)	[7440-31-5]	1000 - 5000	µg/L		
	Titanium (Ti)	[7440-32-6]	80 - 300	µg/L		