

## **Proficiency Testing - Demand, waste water**

|                 |                                   |
|-----------------|-----------------------------------|
| Art. ID         | ERA-578                           |
| Unit            | each                              |
| Deliverydetails | No Dangerous Good /not restricted |

### Description

One 15 mL screw-cap vial yields up to 2 liters after dilution /// Methods: EPA 405.1, EPA 410.1, EPA 410.2, EPA 410.3, EPA 410.4, EPA 415.1, EPA 415.2, EPA 9060, SM 18/19thED 5310B, SM 18/19thED 5310C, SM 18/19thED 5310D, SM 18th ED 5210 B, SM 18th ED 5520 B

| Text/Information | Analyte/Parameter                         | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|---|------------|---------------------|------|--------|--------|
|                  | 5-day BOD Biochemical Oxygen Demand (BOD) |            | 15 - 250            | mg/L |        |        |
|                  | Carbonous BOD (CBOD)                      |            | 15 - 250            | mg/L |        |        |
|                  | Chemical Oxygen Demand (COD)              |            | 30 - 250            | mg/L |        |        |
|                  | Total Organic Carbon (TOC)                |            | 6 - 100             | mg/L |        |        |