

Proficiency Testing - Complex Nutrients, waste water

Art. ID ERA-579
Unit each
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1760 III E1

Description

One 15 mL screw-cap vial yields up to 2 liters after dilution.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Total Kjeldahl-Nitrogen as N		3 - 35	mg/L		
	Total Phosphorus as P		0,5 - 10	mg/L		