

Proficiency Testing - Complex Nutrients, waste water

Art. ID ERA-579
Unit each
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1760 III E1

Description

One 15 mL screw-cap vial yields up to 2 liters after dilution.

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|---------------------------------|------------|---------------------|------|--------|--------|
| | Total Kjeldahl-Nitrogen as N | | 3 - 35 | mg/L | | |
| | Total Phosphorus as P | | 0,5 - 10 | mg/L | | |