

Proficiency Testing - Hardness, waste water

Art. ID ERA-580
Unit each
Deliverydetails No Dangerous Good /not restricted

Description

One 500 mL whole-volume bottle is ready to analyze.

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|---------------------------------|-------------|---------------------|------|--------|--------|
| | Calcium (Ca) | [7440-70-2] | 3,5 - 110 | mg/L | | |
| | Calcium hardness as CaC O3 | | 8,7 - 275 | mg/L | | |
| | Magnesium (Mg) | [7439-95-4] | 4,0 - 40,0 | mg/L | | |
| | Total suspended solids (TSS) | | 23 - 100 | mg/L | | |
| | Total hardness as CaCO3 | | 17 - 440 | mg/L | | |