

## **Proficiency Testing - Minerals, waste water**

Art. ID	ERA-581
Unit	each
Deliverydetails	No Dangerous Good /not restricted

### Description

One 500 mL whole-volume bottle is ready to analyze.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Alkalinity, total (CaCO <sub>3</sub> )		10 - 120	mg/L		
	Chloride (Cl)		35 - 275	mg/L		
	Fluoride (F)		0,3 - 4	mg/L		
	Potassium (K)	[7440-09-7]	4,0 - 40,0	mg/L		
	Sodium (Na)	[7440-23-5]	6 - 100	mg/L		
	Specific Conductance at 25°C		200 - 930	µmhos/cm		
	Sulfate		5 - 125	mg/L		
	Total dissolved solids at 180 °C		140 - 650	mg/L		
	Total solids at 105 °C		140 - 675	mg/L		