

Proficiency Testing - Oil & Grease, waste water

Art. ID	ERA-582
Unit	1 L
Deliverydetails	No Dangerous Good /not restricted

Description

One liter whole-volume glass bottle with a 33-430 thread is ready to analyze. For SPE compatible QC (33-400 thread), use Cat #818 1 liter Boston Round Oil & Grease. Please specify if for SPE compatible PT // Methods: EPA 1664, EPA 413.1, EPA 413.2, EPA 9070, SM 18th ED 5520 B

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Oil and Grease		20 - 100	mg/L		