

## **Proficiency Testing - Total Cyanide, waste water**

Art. ID                      ERA-588  
Unit                         each  
Deliverydetails            No Dangerous Good /not restricted

### Description

One 15 mL screw-cap vial yields up to 2 liters after dilution.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Total cyanide		0,1 - 1	mg/L		
	Total Phenolics (4AAP)		0,5-5	mg/L		
	Amenable Cyanide		0,1 - 1	mg/L		