

## **Proficiency Testing - Total Phenolics (4-AAP), waste water**

Art. ID	ERA-589
Unit	each
Deliverydetails	ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1230 (6.1) II E2

### Description

One 2 mL flame-sealed ampule yields up to 2 liters after dilution. Analyze for total phenolic compounds by 4-AAP methods.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Total phenolics		0,5 - 5	mg/L		
	Phenol	[108-95-2]	0,2 - 5	mg/L		
	2,4-Dinitrophenol	[51-28-5]	0,196 - 5	mg/L		
	2,4-Dichlorophenol	[120-83-2]	0,173 - 5	mg/L		
	2-Chlorophenol	[95-57-8]	0,137 - 5	mg/L		