

## **Proficiency Testing - Halomethanes (THMs), potable water**

Art. ID	ERA-842
Unit	each
Deliverydetails	No Dangerous Good /not restricted

### Description

One 2 mL flame-sealed ampule yields in excess of 200 mL after dilution. Contains all of the analytes below at 5-50 µg/L /// Methods: EPA 502.1, EPA 502.2, EPA 503.1, EPA 524.1, EPA 524.2, EPA 551

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Bromodichloromethane	[75-27-4]	5,0 - 50	µg/L		
	Bromoform (Tribromomethane)	[75-25-2]	5,0 - 50	µg/L		
	Chlorodibromomethane	[124-48-1]	5,0 - 50	µg/L		
	Chloroform (Trichloromethane)	[67-66-3]	5,0 - 50	µg/L		