

## **Proficiency Testing - Chlordane, potable water**

Art. ID ERA-845  
Unit each  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

### Description

One 2 mL flame-sealed ampule yields up to 2 liters after dilution. Each standard contains technical chlordane at 2-20 µg/L.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Chlordane (total)		2,0 - 20,0	µg/L		