

## **Proficiency Testing - EDB/DBCP/TCP, potable water**

Art. ID                      ERA-847  
Unit                         each  
Deliverydetails            No Dangerous Good /not restricted

### Description

One 2 mL flame-sealed ampule yields in excess of 200 mL after dilution. Each lot contains all analytes below at 0.05-2 µg/L /// Methods: EPA 504.1, EPA 551

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	1,2,3-Trichloropropane	[96-18-4]	0,1 - 2	µg/L		
	1,2-Dibromo-3-chloropropane	[96-12-8]	0,1 - 2	µg/L		
	Ethylene Dibromide	[106-93-4]	0,1 - 2	µg/L		