

Proficiency Testing - Dioxin, potable water

Art. ID	ERA-857
Unit	each
Deliverydetails	ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

One 2 mL flame-sealed ampule yields up to 2 liters after dilution. Each standard contains 2, 3, 7, 8-TCDD at 20-100 pg/L.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,3,7,8-Tetrachlorodibe nzo-p-dioxin	[1746-01-6]	20 - 100	pg/L		