

## **Proficiency Testing - Low-Level Total Residual Chlorine**

Art. ID                      ERA-881  
Unit                         each  
Deliverydetails            No Dangerous Good /not restricted

### Description

Designed for testing at low µg/L levels. One 2 mL flame-sealed ampule yields up to 2 liters after dilution. Use with sensitive titrimetric or colorimetric methods.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Total residual chlorine		20 - 250	µg/L		