

## **Proficiency Testing - Sulfide, waste water**

Art. ID                      ERA-891  
Unit                         each  
Deliverydetails            No Dangerous Good /not restricted

### Description

One 10 mL flame-sealed ampule yields up to 1 liter after dilution. Preserved sample is guaranteed stable. Analyze for sulfide by titrimetric or colorimetric methods or ISE.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Sulfide		2,0 - 10	mg/L		