

QuiKResponse - PCBs on Polyurethane foam

Art. ID	ERA-1112QR
Unit	each
Deliverydetails	No Dangerous Good /not restricted

Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 2 mL flame-sealed ampule plus one polyurethane foam (PUF). Contains one Aroclor, randomly selected from the list below, at 1?15 µg/sample after preparation.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aroclor 1016	[12674-11-2]	1,0 - 15	µg/sample		
	Aroclor 1221	[11104-28-2]	1,0 - 15	µg/sample		
	Aroclor 1232	[11141-16-5]	1,0 - 15	µg/sample		
	Aroclor 1242	[53469-21-9]	1,0 - 15	µg/sample		
	Aroclor 1248	[12672-29-6]	1,0 - 15	µg/sample		
	Aroclor 1254	[11097-69-1]	1,0 - 15	µg/sample		
	Aroclor 1260	[11096-82-5]	1,0 - 15	µg/sample		