

## QuiKResponse - Metals on Filter Paper

|                 |                                   |
|-----------------|-----------------------------------|
| Art. ID         | ERA-1125QR                        |
| Unit            | each                              |
| Deliverydetails | No Dangerous Good /not restricted |

### Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One filter paper sample is ready to be analyzed.

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit      | Method | Source |
|------------------|-------------------|-------------|---------------------|-----------|--------|--------|
|                  | Antimony (Sb)     | [7440-36-0] | 30 - 1200           | µg/filter |        |        |
|                  | Arsenic (As)      | [7440-38-2] | 30 - 1200           | µg/filter |        |        |
|                  | Barium (Ba)       | [7440-39-3] | 30 - 1200           | µg/filter |        |        |
|                  | Beryllium (Be)    | [7440-41-7] | 30 - 1200           | µg/filter |        |        |
|                  | Cadmium (Cd)      | [7440-43-9] | 30 - 1200           | µg/filter |        |        |
|                  | Chromium (Cr)     | [7440-47-3] | 30 - 1200           | µg/filter |        |        |
|                  | Cobalt (Co)       | [7440-48-4] | 30 - 1200           | µg/filter |        |        |
|                  | Copper (Cu)       | [7440-50-8] | 30 - 1200           | µg/filter |        |        |
|                  | Lead (Pb)         | [7439-92-1] | 30 - 1200           | µg/filter |        |        |
|                  | Manganese (Mn)    | [7439-96-5] | 30 - 1200           | µg/filter |        |        |
|                  | Nickel (Ni)       | [7440-02-0] | 30 - 1200           | µg/filter |        |        |
|                  | Phosphorous       |             | 30 - 1200           | µg/filter |        |        |
|                  | Selenium (Se)     | [7782-49-2] | 30 - 1200           | µg/filter |        |        |
|                  | Silver (Ag)       | [7440-22-4] | 30 - 1200           | µg/filter |        |        |
|                  | Thallium (Tl)     | [7440-28-0] | 30 - 1200           | µg/filter |        |        |
|                  | Zinc (Zn)         | [7440-66-6] | 30 - 1200           | µg/filter |        |        |