

## **QuiKResponse - Mercury in Impinger Solution**

Art. ID                      ERA-1128QR  
Unit                         each  
Deliverydetails            ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

### Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One impinger solution sample is ready to be analyzed.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Mercury (Hg)	[7439-97-6]	1,0 - 30,0	ng/mL		