

## **QuiKResponse - Hexavalent Chromium in Impinger Solution**

Art. ID ERA-1132QR  
Unit each  
Deliverydetails No Dangerous Good /not restricted

### Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One impinger solution sample is ready to be analyzed.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Hexavalent chromium	[18540-29-9]	50 - 800	µg/L		