

QuiKResponse - Ammonia in Impinger Solution

Art. ID ERA-1145QR
Unit each
Deliverydetails No Dangerous Good /not restricted

Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One impinger solution sample is ready to be analyzed.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Ammonium (NH ₄ ⁺)	[14798-03-9]	0,1 - 10	mg/L		