

## **QuiKResponse - Glycols in Water**

Art. ID	ERA-401QR
Unit	each
Deliverydetails	No Dangerous Good /not restricted

### Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 2 mL flame sealed ampule yields up to 2 L after dilution. Use with EPA methods 8015B, 8430 and 1671.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Diethylene glycol	[111-46-6]				
	Ethylene Glycol	[107-21-1]				
	Propylene glycol	[57-55-6]				
	Tetraethylene glycol	[112-60-7]				
	Triethylene glycol	[112-27-6]				