

QuiKResponse - Solids, Concentrate, waste water (WasteWatR TM)

Art. ID	ERA-4032QR
Unit	each
Deliverydetails	No Dangerous Good /not restricted

Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 23 mL screw-cap vial with a powder yields 1 liter after dilution.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Total suspended solids (TSS)		23 - 100	mg/L		
	Total dissolved solids at 180 °C		140 - 650	mg/L		
	Total solids at 105 °C		140 - 675	mg/L		