

QuiKResponse - Oil & Grease Concentrate, waste water (WasteWatR TM)

Art. ID ERA-4122QR
Unit each
Deliverydetails No Dangerous Good /not restricted

Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 23 mL screw-cap vial yields up to 2 liters after dilution.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Oil and Grease		20 - 100	mg/bottle		