

## **QuiKResponse - Solids, waste water (WasteWatR TM)**

Art. ID	ERA-499QR
Unit	each
Deliverydetails	No Dangerous Good /not restricted

### Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 500 mL whole-volume bottle is ready to analyze.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Total suspended solids (TSS)		20 - 100	mg/L		
	Total Dissolved Solids at 180 °C (TDS)		140 - 800	mg/L		
	Total solids at 105 °C		140 - 800	mg/L		