

## **QuiKResponse - Simple Nutrients, waste water (WasteWatR TM)**

Art. ID                      ERA-505QR  
Unit                         each  
Deliverydetails            No Dangerous Good /not restricted

### Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 15 mL screw-cap vial yields up to 2 liters after dilution.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Ammonia (as N)		0,65 - 19	mg/L		
	Nitrate + Nitrite as N		0,25 - 40	mg/L		
	Nitrate (NO <sub>3</sub> -) as N		0,25 - 40	mg/L		
	Ortho-Phosphate as P		0,5 - 5,5	mg/L		