

QuiKResponse - Simple Nutrients, waste water (WasteWatR TM)

Art. ID ERA-505QR
Unit each
Deliverydetails No Dangerous Good /not restricted

Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 15 mL screw-cap vial yields up to 2 liters after dilution.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Ammonia (as N)		0,65 - 19	mg/L		
	Nitrate + Nitrite as N		0,25 - 40	mg/L		
	Nitrate (NO ₃ -) as N		0,25 - 40	mg/L		
	Ortho-Phosphate as P		0,5 - 5,5	mg/L		