

## **QuiKResponse - Minerals, waste water (WasteWatR TM)**

Art. ID	ERA-506QR
Unit	each
Deliverydetails	No Dangerous Good /not restricted

### Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 500 mL whole-volume bottle is ready to analyze.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Alkalinity (as CaCO <sub>3</sub> )		25 - 400	mg/L		
	Chloride (Cl)		35 - 275	mg/L		
	Fluoride (F)		0,4 - 4	mg/L		
	Potassium (K)	[7440-09-7]	4,0 - 40,0	mg/L		
	Sodium (Na)	[7440-23-5]	10 - 100	mg/L		
	Specific Conductance at 25°C		200 - 1200	µmhos/cm		
	Sulfate		5 - 125	mg/L		
	Total dissolved solids at 180 °C		140 - 800	mg/L		
	Total solids at 105 °C		140 - 800	mg/L		