

QuiKResponse - Hardness, waste water (WasteWatR TM)

Art. ID ERA-507QR
Unit each
Deliverydetails No Dangerous Good /not restricted

Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 500 mL Whole-Volume bottle is ready to analyze.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Calcium (Ca)	[7440-70-2]	3,5 - 110	mg/L		
	Calcium hardness as CaC O3		8,7 - 275	mg/L		
	Magnesium (Mg)	[7439-95-4]	2,0 -40,0	mg/L		
	Total suspended solids (TSS)		23 - 100	mg/L		
	Total hardness as CaCO3		17 - 440	mg/L		