

### **QuiKResponse - Mercury, waste water (WasteWatR TM)**

Art. ID ERA-514QR  
Unit each  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1760 III E1

#### Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 15 mL screw-cap vial yields up to 1 liter after dilution. Analyze for total mercury.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Mercury (Hg)	[7439-97-6]	2,0 - 30,0	µg/L		