

## **QuiKResponse - Solids, Concentrate, potable water (PotableWatR TM)**

Art. ID	ERA-5152QR
Unit	each
Deliverydetails	No Dangerous Good /not restricted

### Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 23 mL screw-cap vial with a powder yields 1 liter after dilution.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Total suspended solids (TSS)		23 - 100	mg/L		
	Total dissolved solids at 180 °C		200 - 450	mg/L		
	Total solids at 105 °C		223 - 550	mg/L		