

QuiKResponse - Demand, waste water (WasteWatR TM)

Art. ID	ERA-516QR
Unit	each
Deliverydetails	No Dangerous Good /not restricted

Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 15 mL screw-cap vial yields up to 2 liters after dilution.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Biochemical Oxygen Demand (BOD)		15 - 250	mg/L		
	Carbonous BOD (CBOD)		15 - 250	mg/L		
	Chemical Oxygen Demand (COD)		30 - 250	mg/L		
	Total Organic Carbon (TOC)		6 - 100	mg/L		