

QuiKResponse - Demand, waste water (WasteWatR TM)

| | |
|-----------------|-----------------------------------|
| Art. ID | ERA-516QR |
| Unit | each |
| Deliverydetails | No Dangerous Good /not restricted |

Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 15 mL screw-cap vial yields up to 2 liters after dilution.

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|---------------------------------|------------|---------------------|------|--------|--------|
| | Biochemical Oxygen Demand (BOD) | | 15 - 250 | mg/L | | |
| | Carbonous BOD (CBOD) | | 15 - 250 | mg/L | | |
| | Chemical Oxygen Demand (COD) | | 30 - 250 | mg/L | | |
| | Total Organic Carbon (TOC) | | 6 - 100 | mg/L | | |