

## **QuiKResponse - Tin & Titanium, waste water (WasteWatR TM)**

Art. ID ERA-517QR  
Unit each  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1789 III E1

### Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 15 mL screw-cap vial yields up to 1 liter after dilution. Use with AA, ICP-OES or ICP-MS methods.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Tin (Sn)	[7440-31-5]	1000 - 5000	µg/L		
	Titanium (Ti)	[7440-32-6]	80 - 300	µg/L		